

Instructions for Use

1. **Turn** patient as far as possible to the side position.



2. **Angle** the cushion so the "tips" point toward bed.




3. **Press Top** upper tip to fit under the scapula (shoulder blade).



4. **Press Bottom** after placing cushion as desired, and press the lower tip under buttocks.



5. **Adjust** to avoid sacrum/coccyx or wound area. See "Hand Check" suggestion below to achieve good support without touching wounds.



6. **Roll Back** on cushion. Best to hold cushion in place during roll.



7. **Pressure Check:** Slide a hand under the cushion to assure proper placement away from pressure points **and** wounds)



Comfort Tips

Every patient is different. To find a position that can achieve maximum comfort and pressure relief, you may alternate the cushion to fit. Place above or below the wound; set closer to the body, or further away; angle the upper tip closer than lower tip or vice versa.

Nursing Advice

We at Jewell Nursing Solutions understand it's not always easy to turn and support a person off the pressure injury properly. If you are having difficulty with turning and repositioning, please see our "Learn to Turn" tutorial for more information. Or feel free to contact us directly if you need specific nursing advice or want to learn about a few tricks of the trade!

www.JewellNursingSolutions.com

(650) 294-8557