



Cavilon™

Skin Care Solutions

3M™ Cavilon™ Advanced Skin Protectant

Protect at-risk periwound skin when using negative pressure wound therapy.

Delivering periwound protection and preventive care.

Periwound skin is at risk of irritants and can delay wound healing, increase discomfort for patients, and have an adverse effect on outcomes. Without the proper management of the periwound skin, negative pressure wound therapy (NPWT) could be compromised resulting in prolonged healing and higher cost of care.



Cavilon Advanced Skin Protectant:



Adheres to wet, weepy, damaged skin creating an environment that allows healing.¹



Does not diminish adherence of 3M acrylic adhesive-based tapes and drapes.



Is breathable, so moisture can evaporate.



Has a low dermatitis potential.



Take a proactive approach to maintaining, managing and protecting at-risk skin under NPWT.

Precautions: The product can increase the adhesion of some adhesive products (e.g., 3M acrylic negative pressure drapes), particularly in the first few days of use. When removing an adhesive product, it is important to exercise care. Use an adhesive remover containing HMDS to remove an adhesive product (including 3M acrylic negative pressure drapes) that has been placed over the film. This is especially important for patients with fragile skin.

Unique Solutions for Skin Protection.

3M™ Cavilon™
Advanced Skin Protectant



3M™ Cavilon™
No Sting Barrier Film



Recommended for use on at-risk or broken, wet/weepy, or denuded periwound skin under NPWT dressings.



Product is housed in sterile packaging.



Use for routine protection of periwound under NPWT dressings.



Use for routine protection of periwound skin.



*Applies to wipes and swabs only

3M™ Cavilon™ Advanced Skin Protectant is available in two sizes to suit a variety of coverage needs.



3M™ Cavilon™ Advanced Skin Protectant 5050



3M™ Cavilon™ Advanced Skin Protectant 5051



Ordering Information

Cat. No	Size	Approx. Coverage	Ideal for	Items/Box
5050	2.7ml applicator	100 sq. in.	<ul style="list-style-type: none"> Periwound skin management under NPWT applications 	20
5051	0.7ml applicator	25 sq. in.	<ul style="list-style-type: none"> Peristomal, peri-fistula, peri-tube and periwound skin damage protection and management Protect skin around IV insertion site 	20

For more information on Cavilon Advanced Skin Protectant, go to [3M.com/Cavilon](https://www.3M.com/Cavilon) or contact your 3M account representative.



3M Company
2510 Conway Ave
St. Paul, MN 55144 USA

Phone 1-800-275-4524 (NPWT products)
1-800-228-3957
Web 3M.com/medical

Note: Specific indications, limitations, contraindications, warnings, precautions and safety information exist for these products and therapies some of which may be Rx only. Please consult a clinician and product instructions for use prior to application.

1. Acton C, Ivins N, Bainbridge P, Browning P. Management of incontinence-associated dermatitis patients using a skin protectant in acute care: a case series. J Wound Care. 2020;29(1):18-26. doi:10.12968/jowc.2020.29.1.18

© 2023 3M. All rights reserved. 3M and the other marks shown are marks and/or registered marks. Unauthorized use prohibited. Under license in Canada. GL_70-2013-1571-3