


# Cutimed® Sorbion®

## A superabsorbent dressing featuring multiple features and benefits

**New Product Sizes Available!†**

Cutimed® Sorbion® superabsorbent dressings, featuring **Hydration Response Technology™**, help promote the healing needs of chronic wounds using a unique combination of high absorption, retention capacity, and soft debridement.



### Water-retaining polymers and fiber matrix

- Wound exudate is safely absorbed and locked in vertically

### Special three-dimensional wound contact layer


- Minimizes the contact area with the wound bed for an atraumatic dressing change

### Special pore layer

- Supports the transfer of wound exudate to the inner core


### Superabsorber

- Reduction in risk of potential skin maceration
- For moderate to high exudate levels



†New product sizes for Cutimed® Sorbion® Sachet will be available on September 1, 2024.


## Offering a variety of shapes and sizes\*



### Cutimed® Sorbion® Sachet

Easy application – either side of the dressing may be placed on wound

Size*	Order Code
3" x 3"	7323200
4" x 4"	7323206
<b>NEW SIZE</b> 5" x 5"†	<b>7323268†</b>
4" x 8"	7323209
<b>NEW SIZE</b> 6" x 8"†	<b>7323269†</b>
8" x 8"	7323215
12" x 8"	7323218



### Cutimed® Sorbion® Border

Includes a high tack silicone border for atraumatic removal

Size*	Order Code
4" x 4"	7101000
6" x 6"	7101002
8" x 8"	7101003
Sacrum 8" x 9"	7101008

\*All items sold in boxes of 10 dressings.

**Reimagine what your wound dressings can do.**

Connect with your local Essity sales representative to learn more.

BSN Medical Inc., an Essity company  
 5825 Carnegie Blvd. Charlotte, NC 28209-4633  
 Tel. (+1) 704 554 9933 Fax: (+1) 800 835 4325  
 To order toll-free: (+1) 800 537 1063