NOTICE:

Effective Immediately



SilvaKollagen® Gel is a silver antimicrobial gel with hydrolyzed collagen that begins working immediately on contact. The patented formula absorbs 8 to 10 times its weight in fluid and contributes moisture if the wound is dry. Published abstracts report improved healing^{1,4}, improved pain control with decreased need for narcotics^{1,4} and expedient debridement³. SilvaKollagen® Gel also helps control bleeding, deodorizes, and minimizes scarring.

Learn more at DermaRite.com or call 800-337-6296.



- Publications:

 1. Heasley D. Evaluation of a silver collagen-based amorphous gel for the use of treatment of a deep partial-thickness burn. *J WOCN.* 2009;36(3S):S9.

 2. Khandaker K, Kohl DS. Evaluation of the use of a silver collagen amorphous gel in the healing of postsurgical and dehisced lesions. *J WOCN.* 2009;36(3S):S13.

 3. Heasley D. Evaluation of a silver collagen based amorphous gel in expedience of wound debridement. *J WOCN.* 20011;38(3S):S18.

 4. Heasley D. Evaluation of a silver collagen based gel and attained pain relief in arterial ulcers. *J WOCN.* 20011;38(3S):S18.